


## bio-m® AGITATORS

### GENERAL INFORMATION ON AGITATOR SIZING



- PROCESS LAY-OUT
- MECHANICAL DESIGN



BIO- AND PROCESSTECHNOLOGY | STERILE PLANT CONSTRUCTION | AUTOMATION



## INTRODUCTION:



There are calculations available for sizing an agitator, with regard to the process layout as well as the mechanical features. The calculations are similar for magnetically coupled and mechanically sealed agitators.

On the following pages, zeta has put together basic information on these topics:

Agitator Sizing	Page 3
Viscosity	Page 4
Agitation Task: Homogenisation	Page 5
Agitation Task: Suspending	Page 6
Agitation Task: Heat Transfer	Page 7
Agitation Task: Gas Dispersion	Page 8
Typical Impellers	Page 9
Calculation of the required Motor Power	Page 10
Calculation of the necessary shaft diameter	Page 11
First critical speed	Page 12
Questionnaire	Page 13

If there are any questions don't hesitate to contact us



## AGITATOR SIZING

### 1. BASIC INFORMATION

about the process, vessel dimensions, operating conditions etc.

### 2. INFORMATION ABOUT THE PRODUCT

at least: S.G., Viscosity (refer to data sheet)

### 3. AGITATION TASK

- Blending
- Suspending of Solids
- Heat Transfer
- Gas Dispersion etc.

### 4. SIZING TO THE PROCESS REQUIREMENTS

(by experience, CFD, lab-scale tests, scale up)

Result: Definition of impeller design, - diameter, -speed

### 5. CALCULATION OF THE POWER REQUIREMENT / TORQUE

Equations:  $P \sim s.g. \times n^3 \times d^5$ ,  $Mt = P / (2 \times \pi \times n)$

### 6. MECHANICAL DESIGN

i.e. type of shaft bearing, shaft diameter, critical speed, lantern, type and size of drive system, magnetic coupling etc.

### 7. FURTHER CHARACTERISTICS

for instance size of vessel flange, m.o.c., surface finish etc.



## VISCOSITY

Viscosity is a characteristic measure for the internal friction of a fluid.

Usually the dynamic viscosity  $\eta$  is given, which correlates to the kinetic viscosity  $\nu$  and the specific gravity  $\rho$  :  $\eta = \rho \times \nu$

Standard dimension of the dynamic viscosity: mPa s  
 (milli Pascal Sekond) = 1mPa s= 1Pa s/1000, 1 Pa s = 1 N s / m<sup>2</sup>

Old dimension (should not be used any longer):  
 cP (centi Poise), 1 cP = 1 mPa s

The viscosity strongly depends on the temperature; a viscosity therefore should always be given together with the temperature, at which the measurement was done.

With **Newtonian** fluids, the viscosity depends on the temperature only. With **non-Newtonian** fluids there are more influences, especially shear forces (shear thinning) and/or duration of the shear stress (thixotrope)

### EXAMPLES:

PRODUCT	TEMPERATURE	VISCOSITY
	°C	mPa s
Solvents (Aceton, Toluol etc.)	20	< 1
Water	20	≈ 1
Liqueur	20	12
Sugar solution 50%	20	15
Salad oil	20	85
Sugar solution 64%	20	120
Egg liqueur	20	620
Glycerin 100%	30	624
Glycerin 100%	25	945
Glycerin 100%	20	1500
Yoghurt	20	≈ 1000
Honey	45	≈ 2000
Apple puree	20	≈ 10000



## AGITATING TASK: BLENDING / HOMOGENISATION

### 1. DEFINITION

Blending of miscible liquids, within a certain, given or required mixing time, to a certain, given or required degree/quality.

### 2. MIXING TIME DIAGRAM

This diagram is produced in a lab- or pilot scale test procedure and is valid for a certain impeller type and mounting situation - in any vessel size provided there is geometric similarity (impeller shape, d/D ratio, H/D-ratio etc.)

In this test procedure the mixing time  $\theta$  is measured which it takes to achieve a certain mixing quality (usually 95 % of pH value, of temperature etc.), with different speeds of rotation  $n$  of the impeller.

The diagram shows the product  $n \times \theta$  over Reynolds No.  $Re$

### 3. MEASURING METHODS

There are different methods to determine the mixing time which may deliver slightly different results.

In the lab scale, the de-colouring method is used frequently: A blue iodine-starch solution is added to the liquid. Then a sodiumthiosulfate solution is added and the time is measured how long it takes the agitator to decolour the liquid.

Alternative methods can be to measure the time it takes to achieve an even conductivity or temperature.

### 4. EXAMPLE

Marine type propeller, speed  $n = 200$  1/min

$Re > 10\,000$  (turbulent flow)

$n \times \theta = 100$  (from Mixing time diagram)

Therefore:

$\theta = (100 / 200)$  min, mixing time is 30 sec

(provided there is geometric similarity to the test scale)



## AGITATING TASK: SUSPENDING OF SOLIDS

1. CALCULATION OF THE SETTLING RATE OF A PARTICLE (SIZE  $d_p$ ):

$$v_s = f(\Delta\rho, \eta, d_p)$$

(with the help of the Reynolds-, Archimedes- und Nusselt-Number)

2. CALCULATION OF THE SETTLING RATE OF A PARTICLE SWARM:

$$v_{ss} = v_s (1 - c_v)^x \quad (c_v = \text{solids content})$$

3. CALCULATION OF THE SPECIFIC SETTLING POWER:

$$E_{ss} = (\Delta\rho / \rho_{\text{liquid}}) \times c_v \times g \times v_{ss}$$

4. CALCULATION OF THE NECESSARY AGITATION POWER P:

- $P > E_{ss}$
- $P/V = f$  (tank diameter  $d_1$ , expected suspension quality)  
( $P/V$  = specific power input, power input per  $m^3$ )
- Interrelationship between  $P/V$  and  $d_1$  is not linear;  
(The curve  $P/V$  over  $d_1$  has been determined by Zehner, Kipke et al.)



## AGITATING TASK: HEAT TRANSFER

### 1. FLOW OF HEAT Q:

Q is proportional to the heat transfer area A and to the difference in temperature:

$$Q = k \times A \times \Delta T$$

### 2. HEAT TRANSFER COEFFICIENT K:

$$1/k = 1/\alpha_i + s/\lambda_w + 1/\alpha_a$$

$\alpha_i$  = Internal heat transfer coefficient (Product -> Vessel wall)

s = Thickness of vessel wall [m]

$\lambda_w$  = Heat conductivity of the vessel wall [W / mK]

$\alpha_a$  = External heat transfer coefficient [W / m<sup>2</sup>K]

### 3. CALCULATION OF $\alpha_i$ :

3.1 Prandtl-No. (viscosity, heat capacity and -conductivity of the product) :

$$Pr = \eta_R \times c_p / \lambda_R$$

3.2 Nusselt-No.

$$Nu = C \times Re^{2/3} \times Pr^{1/3}$$

(if Re >200 ; neglecting the difference in viscosity near the tank wall)

C = f ( typ of impeller, mounting)

C = 0.4 .....0.9

3.3 Internal heat transfer coefficient

$$\alpha_i = Nu \times \lambda_R / d_1, \quad \alpha_i \sim P^{2/9}$$

(d<sub>1</sub> = Tank Diameter, P = agitator power)



## AGITATING TASK: GASSING

### 1. PURPOSE OF MIXING

to intensify mass transfer by means of

- Magnification of the specific mass transfer surface  $a$  /  
minimization of the bubble size
- max. concentration grade  $\Delta c$  by high turbulence
- Magnification of gas hold up time
- short mixing time

### 2. CHARACTERISTIC FIGURES

gassing rate: gas input (volume) per minute / vessel  
working volume

$$vvm = q / V [1/\text{min}]$$

Gas velocity: gas input (volume) per minute /  
cross sectional area of the tank

$$v_L = q \times 4 / d_1^2 \times \pi$$

specific gas throughput:  $Q = q / n_2 \times d_2^3$

Froude-No.:  $Fr = n_2^2 \times d_2 / g$

### 3. CALCULATION

the calculation has to be done by the iteration method:

1. Specification of impeller-type, -quantity, -diameter
2. Specification of the impeller speed of rotation (by  
experience or looking for a certain P/V value)
3. Determination of the Power No. „Ne“ from the diagram  
„Ne over Q“, Froud No. as parameter
4. Check: the impeller must not be flooded by the amount of  
gas (diagram, specific for each impeller type)
5. Depending on 4., a new calculation may be necessary,  
taking a different impeller diameter, different speed;  
Optimization of the process lay out
6. Calculation of the power requirement and spec. power P/V
7. Calculation of the mass transfer:

$$k_L a = C \times (P/V)^A \times (v_L)^B$$

coefficients

$A = 0,3 \dots 0,8$  (depending on the viscosity)








$B = 0,5 \dots 0,3$  (coalescenc .... no coalescence )

proportional factor

$C = f$  ( impeller ) (e.g. Rushton turbine:  $C = 0.41$ )



## TYPICAL IMPELLERS:

TYPE	TYPICAL AGITATION TASK, CHARACTERISTICS	VISCOSITY RANGE	DIAMETER RATIO	
		MPaS	D/D	
Marine type propeller	Homogenise Suspend solids circulate; Axial flow pattern, High rpm	1 - 1000	0,15 - 0,35	
Pitched blade Impeller	Homogenise Suspend solids Improve heat transfer; Axial flow pattern, with a radial component; medium rpm	1 - 2500	0,2 - 0,4	
Rushton Turbine	Gassing; Radiale flow pattern; Usually multi stage; High speed (rpm) / high power input	1 - 1000	0,3 - 0,45	
Segment Impeller	Gentle homogenisation, improve heat transfer; Axial flow pattern, with a radial component; Usually multi stage Low speed	1 - 1000	0,4 - 0,5	
Cross beam Impeller	Homogenise Suspend solids Improve heat transfer; Axial flow pattern; Usually multi stage	100 - 50000	0,4 - 0,8	
Counterflow Impeller	Homogenise, gassing, Suspend solids, Improve heat transfer; Axial counter flow (inside down, outside up); Always multi stage	100 - 50000	0,5 - 0,9	
U/Z - Impeller	Homogenise Improve heat transfer; Always multi stage; Often used in food industry	100 - 10000	0,4 - 0,6	

You may also refer to our brochure and to DIN 28131



## CALCULATION OF THE REQUIRED MOTOR POWER

### 1. PRODUCT CHARACTERISTICS

at least: specific gravity  $\rho$  [ kg/m<sup>3</sup>], viscosity  $\eta = f ( T )$  [ Pa s]

### 2. DESIGN FEATURES ACCORDING TO THE PROCESS:

Type and number of impellers

Impeller diameter  $d_2$  [m]

Rotational speed  $n_2$  [1/s]

### 3. CALCULATION OF THE RYNOLDS NUMBER RE:

$$Re = \frac{n_2 \times d_2^2 \times \rho}{\eta}$$

### 4. DETERMINATION OF THE POWER NO. NE (NEWTON-NO.):

$Ne = f ( Re, \text{Impeller Type, diameter ratio etc.})$

to be taken from the relevant diagram Power No. vs Reynolds No.

### 5. CALCULATION OF THE POWER REQUIRED BY THE IMPELLER:

$$P_v = Ne \times \rho \times n_2^3 \times d_2^5$$

### 6. CALCULATION OF THE TORQUE:

$$M_t = \frac{P}{2 \times \pi \times n_2}$$

### 7. FACTOR FOR ANY LOSSES, SAFETY ETC.

$$P_{\text{Motor}} = 1,1 \dots 1,3 \times P$$



## CALCULATION OF THE NECESSARY SHAFT DIAMETER

### 1. DETERMINATION OF THE RADIAL FORCE AT THE IMPELLER (turbulent flow pattern)

$$F_R = c_R \times \rho \times n_2^2 \times d_2^4$$

$n_2$  = speed of rotation,  $d_2$  = impeller diameter  
 $c_R$  is specific for the impeller type, has to be measured,  
 $c_R = f$  ( impeller, mounting position, first critical speed )  
 $\rho$  = density

### 2. CALCULATION OF THE BENDING MOMENT

$$M_b = F_R \times l, \quad l = \text{shaft length}$$

### 3. CHECK: DOES THE IMPELLER RUN DURING FILLING/ EMPTYING OF THE TANK

if yes:  $M_{b, D} > M_b$

### 4. CALCULATION OF THE MOMENT OF COMPARISON (Hypothesis of the maximum shape changing work)

$$M_v = \sqrt{M_b \cdot D^2 + 0.37 \times M_t^2}, \quad M_t = \text{max. drive torque}$$

### 5. CALCULATION OF THE NECESSARY SECTION MODULUS AGAINST BENDING

$$W_b = M_v / \sigma_v, \quad \sigma_v = \text{max. equivalent stress}$$

### 6. CALCULATION OF THE SECTION MODULUS OF A SHAFT

$$W_{b, W} = (\pi \times d_W^3) / 32, \quad W_{b, W} \text{ must be } > W_b$$




## FIRST CRITICAL SPEED

Any agitator is an oscillating system and therefore has a natural frequency or resonance frequency, called the first critical speed.

The agitator's speed of rotation must keep away (at least 25 %) from this frequency. This applies especially for speed controlled drives (VFD).

Otherwise, the shaft will vibrate more and more and finally will break down.

The first critical speed can be determined quite exactly by means of a field transfer matrix which contains all major agitator components:



# QUESTIONNAIRE FOR BIO-M<sup>®</sup> AGITATORS

Inquiry from: \_\_\_\_\_ Inquiry No: \_\_\_\_\_  
 worked out by: \_\_\_\_\_ Tel.: \_\_\_\_\_ Fax.: \_\_\_\_\_

PRODUCT	QUANTITY	S.G.	VISCOSITY	PARTICLE SIZE	TEMP.
	kg	kg/m <sup>3</sup>	mPa s	µm	°C
Components					
FINAL PRODUCT					

PRODUCT CHARACTERISTICS:	
Flow behaviour:	<input type="checkbox"/> newtonian <input type="checkbox"/> non-newtonian
agitation task:	<input type="checkbox"/> homogenize/mix <input type="checkbox"/> suspend solids
	<input type="checkbox"/> heat transfer <input type="checkbox"/> dissolve gas
	<input type="checkbox"/> disperse liquid-liquid
intensity:	
Production process:	<input type="checkbox"/> batchwise <input type="checkbox"/> continuously
	flow rate [m <sup>3</sup> /h]:
VESSEL:	
Please attach a vessel drawing and a nozzle list, if necessary	
	<input type="checkbox"/> already existing <input type="checkbox"/> to be supplied by zeta <input type="checkbox"/> to be supplied by others
Design:	<input type="checkbox"/> cylindrical, vertical <input type="checkbox"/> cylindrical, horizontal <input type="checkbox"/> rectangular <input type="checkbox"/> with isolation
diameter/length [mm]:	height [mm] zyl.:                      total:
kind of bottom:	kind of top:
filling volume [Ltr.] from:	to:
size of manhole:	
	<input type="checkbox"/> agitator runs during filling or emptying
equipment inside the vessel (baffles, heating coils etc.):	
working pressure [bar] from:	to:                      temperatur [°C] from:                      to:
external heat exchanging system:	<input type="checkbox"/> double jacket <input type="checkbox"/> half pipe
	<input type="checkbox"/> others:
agitator support:	DN                      PN
	acc.to: <input type="checkbox"/> nozzle length [mm]:

BIO- AND PROCESSTECHNOLOGY | STERILE PLANT CONSTRUCTION | AUTOMATION



material product wetted:	others	<input type="checkbox"/> surface finish inside:
	outside:	
Ex-Zone [0,1 or 2] inside vessel:	outside vessel:	Temperature class: T
<b>DRIVE:</b>		
purchased by customer	<input type="checkbox"/> yes	<input type="checkbox"/> no
	<input type="checkbox"/> make preferred:	
fixed speed:	<input type="checkbox"/> Motor pole-changeable	<input type="checkbox"/> variable speed
	<input type="checkbox"/> by means of	
tension [V]:	frequency [Hz]:	protection IP
	expl.-proof:	
additional information:		
<b>AGITATOR DESIGN:</b>		
	<input type="checkbox"/> hygienic	<input type="checkbox"/> sterile
	<input type="checkbox"/> CIP	<input type="checkbox"/> SIP
Is a foot steady bearing - to support the shaft - permitted?		
	<input type="checkbox"/> yes	<input type="checkbox"/> no
shaft sealing:	<input type="checkbox"/> not necessary	<input type="checkbox"/> seal ring
	<input type="checkbox"/> mechanical seal	<input type="checkbox"/> magnetic coupling
material of O-rings:	<input type="checkbox"/> FDA conforming	
purchased by customer?	<input type="checkbox"/> yes	<input type="checkbox"/> no
	<input type="checkbox"/> make preferred:	
additional seal equipment required:		
Impeller:	<input type="checkbox"/> zeta's choice	<input type="checkbox"/> customer's experience:
Additional remark:		

Fax: +41/55/460 13 33

