

Rubber Fab Technologies Group's

LASER ENGRAVED HYGIENIC SEALS



Laser Engraved Hygienic Seals

- Laser Engraved Hygienic Seals, Screens and Orifice Plates insure permanent batch and/or lot traceability.
- Engraved compound information can be used to replace color dot compound code.
- Laser engraving available on:

Buna	EPDM
Silicone	FKM
PTFE	Tuf-Steel®
Tuf-Flex®	Compression Control
Orifice Plates	Screens
- Laser engraving available on hygienic seals from 1/2" through 4.0" sizes.
- Engraving will remain readable throughout the seal's service life.
- Special application markings can be engraved on seal at customer's request.
- Consult factory for limitations on etching space and pricing.
- Procedure usually adds one day to ship time, depending on quantity.
- Metal Detectable Elastomers available.
- Animal Derived Ingredient Free.




© Torque-Rite, Tuf-Steel, Tuf-Flex and ADI Free are registered trademarks of Rubber Fab Technologies Group.

Controlled Compression Gasket: United States Patent Number 6,857,638.

Torque-Rite: United States Patent Number 6,082,941.

Rubber Fab Technologies Group's Metal Detectable products are patent pending.

* Buna is not an ADI Free elastomer.



Rubber Fab Technologies Group
Limecrest Road

P.O. Box 626 • Andover, NJ 07821
phone: 973-579-2959 fax: 973-579-7275

www.rubberfab.com




For Additional Information...

For detailed information on any Rubber Fab hygienic seal or product line and ordering information, contact Rubber Fab Technologies Group directly at (973) 579-2959 for your local distributor. For literature, email your requests to sales@rubberfab.com or download a specific brochure at www.rubberfab.com.



In Line System Validation
Smart Gasket
#RF-100



Fluid Conditioning Technology
Screen Hygienic Seals
#RF-110



High Purity
Platinum Silicone Hygienic Seals
#RF-170



Concentric and Eccentric
Orifice Plates
#RF-130



Hygienic Seal
Material Guidelines
#RF-140



The Perfect Surface
Hygienic Seal System
Torque-Rite®
#RF-150

Rubber Fab Technologies Group is a member of:



International Association of
Food Industry Suppliers



Distributed by: